

<b>Name: Purified recombinant protein of Human cytotoxic T-lymphocyte-associated protein 4 (CTLA4), transcript variant 1</b> <b>Product Data Sheet - PROTEIN</b>	<b>Catalog: TP723874</b>
---	--------------------------

**Gene Name:** cytotoxic T-lymphocyte associated protein 4

**GeneBank accession:** NM\_005214

**Description:** Purified recombinant protein of Human cytotoxic T-lymphocyte-associated protein 4 (CTLA4), transcript variant 1 (TP723874)

**Amount:** 10 ug

**Buffer and Storage:** 1 x PBS

**Source:** Recombinant protein was produced with Human CTLA-4, the region of Ala37-Phe162, from gene Accession# NM\_005214

**Protein Accn:** NP\_005205

**Gene Synonym:** ALPS5 | CD | CD152 | CELIAC3 | CTLA-4 | GRD4 | GSE | IDDM12

**Predicted molecular weight:** 39.9 kDa

**Tags:** Fc chimera

**Purity:** >95%, as determined by Coomassie stained SDS-PAGE

**Concentration:** Protein formulated and vialled at >50 ug/mL in specified buffer

\* More data may be available on our website: [TP723874](#)