

Rabbit antibody to c Myc tag

Code	OSM00067W
ID Tag	Rb924-170509-WS
Unit size	150 µl
Immunogen	Synthetic peptide {EQKLISEEDL (410-419)} from human c-Myc (c Myc tag) conjugated to blue carrier protein was used as the antigen.
Conjugate	Unconjugated antibody
Also known	Myc proto-oncogene protein, c-Myc, Transcription factor p64
Host	NZ white rabbit
Purity	Whole serum
Clonality	Polyclonal
Isotype	Polyclonal, whole serum
Applications	IHC, WB using recombinant protein expressing {EQKLISEEDL (410-419)}. It is not recommended for tissue lysates. A dilution of 1 : 300 to 1 : 2000 is recommended. The optimal dilution should be determined by the end user.
Specificity	Specific for c-Myc.
Spcs X-react.	Reacts with c-Myc (EQKLISEEDL).
Format	Lyophilised
Reconstitution	Reconstitute in 150 µl of sterile water. Centrifuge to remove any insoluble material.
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Expiry Date	12 months after reconstitution
Shipping	This item will be shipped to you at ambient temperature in a lyophilised form.
Limitation	For research use only

Related Products

OSM00068W	Sheep antibody to c Myc tag
OSM00075G	Rabbit antibody to c Myc tag