

## Rabbit antibody to c Myc tag

<b>Code</b>	OSM00075G
<b>ID Tag</b>	Rb924-170509-G
<b>Unit size</b>	500 µg
<b>Immunogen</b>	Synthetic peptide {EQKLISEEDL (410-419)} from human c-Myc (c Myc tag) conjugated to blue carrier protein was used as the antigen.
<b>Conjugate</b>	Unconjugated antibody
<b>Also known</b>	Myc proto-oncogene protein, c-Myc, Transcription factor p64
<b>Host</b>	NZ white rabbit
<b>Purity</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	Polyclonal, IgG
<b>Applications</b>	IHC, WB using recombinant protein expressing {EQKLISEEDL (410-419)}. It is not recommended for tissue lysates. A concentration of 10-50 µg/ml is recommended. The optimal concentration should be determined by the end user.
<b>Specificity</b>	Specific for c-Myc.
<b>Spcs X-react.</b>	Reacts with c-Myc (EQKLISEEDL).
<b>Format</b>	Lyophilised
<b>Reconstitution</b>	Reconstitute in 500 µl of sterile water. Centrifuge to remove any insoluble material.
<b>Storage</b>	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
<b>Expiry Date</b>	12 months after reconstitution
<b>Shipping</b>	This item will be shipped to you at ambient temperature in a lyophilised form.
<b>Limitation</b>	For research use only

---

### Related Products

OSM00068W	Sheep antibody to c Myc tag
OSM00067W	Rabbit antibody to c Myc tag