

Rabbit antibody to TRPV1

Code	OSR00246W
ID Tag	Rb254-251107-WS
Unit size	100 µl
Immunogen	A synthetic peptide from the 4th cytoplasmic loop of rat TRPV1 (VR1, capsaicin receptor, vanilloid receptor 1) conjugated to an immunogenic carrier protein was used as the immunogen.
Conjugate	Unconjugated antibody
Also known	Transient receptor potential cation channel subfamily V member 1, TrpV1, osm-9-like TRP channel 1, OTRPC1, Vanilloid receptor 1, Vanilloid receptor type 1-like, Capsaicin receptor, transient receptor potential cation channel, subfamily V, member 1, Vr1, Vr111, VR.5'sv
Host	NZ white rabbit
Purity	Whole serum
Clonality	Polyclonal
Isotype	Polyclonal, whole serum
Applications	IHC, WB. A dilution of 1 : 200 to 1 : 2000 is recommended. The optimal dilution should be determined by the end user.
Specificity	Specific for TRPV1.
Spcs X-react.	Rat. Other species not yet tested.
Format	Lyophilised
Reconstitution	Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material.
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Expiry Date	12 months after reconstitution
Shipping	This item will be shipped to you at ambient temperature in a lyophilised form.
Limitation	For research use only

Related Products

OST00045A Rabbit antibody to TRPV1