

Rabbit antibody to zebrafish TRPV4 (20-70)

Code	OST00071W
ID Tag	Rb437-050608-WS
Unit size	100 µl
Immunogen	A synthetic peptide from aa region 20-70 zebrafish TRPV4, conjugated to an immunogenic carrier protein was used as the antigen.
Conjugate	Unconjugated antibody
Also known	Transient receptor potential cation channel, subfamily V, member 4
Host	NZ white rabbit
Purity	Whole serum
Clonality	Polyclonal
Isotype	Polyclonal, whole serum
Applications	IHC, WB. A dilution of 1: 300 to 1: 2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Specificity	Specific for TRPV4.
Spcs X-react.	Zebrafish. Other species not yet tested.
Format	Lyophilised
Reconstitution	Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material.
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Expiry Date	12 months after reconstitution
Shipping	This item will be shipped to you at ambient temperature in a lyophilised form.
Limitation	For research use only

Related Products

OST00072G	Rabbit antibody to zebrafish TRPV4 (20-70)
OST00417W	Rabbit antibody to zebrafish TRPV4 (800-841)