

Rabbit antibody to TRPM8

Code	OST00144W
ID Tag	Rb544-290608-WS
Unit size	100 µl
Immunogen	A synthetic peptide from the 1 st cytoplasmic loop of chicken TRPM8 (Menthol receptor) conjugated to an immunogenic carrier protein was used as the immunogen.
Conjugate	Unconjugated antibody
Also known	Menthol receptor
Host	NZ white rabbit
Purity	Whole serum
Clonality	Polyclonal
Isotype	Polyclonal, whole serum
Applications	IHC, WB. A dilution of 1 : 300 to 1 : 3000 is recommended. The optimal dilution should be determined by the end user.
Specificity	Specific for TRPM8.
Spcs X-react.	Chicken. Other species not yet tested.
Format	Lyophilised
Reconstitution	Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material.
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Expiry Date	12 months after reconstitution
Shipping	This item will be shipped to you at ambient temperature in a lyophilised form.
Limitation	For research use only