

Rabbit antibody to TRPV1

Code	OST00382W
ID Tag	Rb1868-030912-WS
Unit size	100 µl
Immunogen	A synthetic peptide from aa region 600-660 of chicken TRPV1 (Vanilloid receptor-like protein, VR1, capsaicin receptor) conjugated to an immunogenic carrier protein was used as the antigen.
Conjugate	Unconjugated antibody
Also known	Vanilloid receptor-like protein, VR1, capsaicin receptor, transient receptor potential cation channel, subfamily V, member 1
Host	NZ white rabbit
Purity	Whole serum
Clonality	Polyclonal
Isotype	Polyclonal, whole serum
Applications	IHC, WB (confirmed by recombinant protein). A dilution of 1 : 300 to 1 : 2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Specificity	Specific for TRPV1.
Spcs X-react.	Chicken. Other species not yet tested.
Format	Lyophilised
Reconstitution	Reconstitute in 100 µl of sterile water. Centrifuge to remove any insoluble material.
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Expiry Date	12 months after reconstitution
Shipping	This item will be shipped to you at ambient temperature in a lyophilised form.
Limitation	For research use only

Related Products

OST00105W Rabbit antibody to TRPV1