

Wnt3a alternative peptide (β -catenin pathway agonist)

Code	PG-008
MW (Acetate):	5099.63
Purity:	$\geq 95\%$ (HPLC)
Amount:	10 μg

Formulation: This vial contains 10 μg of the titled compound, which has been lyophilized from 0.2 μm filtered DMSO solution.


Reconstitution: Reconstitute at 20~40 μM in DMSO. Centrifuge the vial before opening.

<Example>

40 μM solution prepared when 10 μg peptide reconstituted by 49.0 μL DMSO.

Storage: Store at -20°C or below for the container vial unopened. Recommend to use the peptide as quickly as possible after opening and reconstitution.

Activity: Measured by its ability to induce TCF/LEF reporter activity



This product is "For Research Use Only"