


StarLight™ Calibration Slides and Kit

Catalog Numbers: 25441, 25442, 25443, 25444, 25445

DESCRIPTION


StarLight™ Calibration Slides feature **vibrant ~6µm fluorescent microspheres** dyed with a single fluorophore for basic imaging checks and calibration. Our four standard versions are appropriate for use with common microsphere filter sets: **Glacial Blue** (360, 450), **Dragon Green** (480, 520), **Envy Green** (525, 565), and **Flash Red** (660, 690), and are available individually or as the full 4-slide collection.




Images of StarLight™ Calibration Slides. The left image is Glacial Blue. The right images are (from top to bottom) Dragon Green, Envy Green, and Flash Red.

SPECTRA INFORMATION


Glacial Blue (360, 450)




Dragon Green (480, 520)



Envy Green (525, 565)



Flash Red (660, 690)



STORAGE AND STABILITY

Storage Store at room temperature. Protect from light and moisture.

Stability StarLight™ slides do not have a formal expiration date. Longevity will depend on proper storage, as well as specific imaging conditions.

This product is for research use only and is not intended for use in humans or for *in vitro* diagnostic use.

ORDERING INFORMATION

Cat. #	Description	Sizes
25441	StarLight™ Calibration Slide – Glacial Blue	1 slide
25442	StarLight™ Calibration Slide – Dragon Green	1 slide
25443	StarLight™ Calibration Slide – Envy Green	1 slide
25444	StarLight™ Calibration Slide – Flash Red	1 slide
25445	StarLight™ Collection – Slide 4-Pack	4 slides

TO ORDER

In The U.S. Call: 1(800) 523-2575 • (215) 343-6484
In The U.S. Fax: 1(800) 343-3291 • (215) 343-0214

In Germany Call: +(49) 06201-845200
In Germany Fax: +(49) 06201-8452020

In Asia Call: (886) 2 8712 0600
In Asia Fax: (886) 2 8712 2677

Order online anytime at www.polysciences.com