

TLQP-21 (Human)

Catalog#	Standard size	Price	Order
003-90	100 µg	\$160.00	<input type="button" value="ORDER"/>

Product description

More information

References

[TLQP and other VGF-derived Peptides](#)

[Orphan Bioactive Peptides from VGF and PC2](#)

[Alzheimer biomarker](#)

Catalog #	003-90
Standard Size	100 µg
Sequence	Thr - Leu - Gln - Pro - Pro - Ser - Ala - Leu - Arg - Arg - Arg - His - Tyr - His - His - Ala - Leu - Pro - Pro - Ser - Arg
Species	Human
M.W	2490.86
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥90%

Thursday 20 December, 2012

11979102 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.