

TLQP-11 / VGF, Prepro (556-566) (Rat)

| Catalog# | Standard size | Price | Order |
|----------|---------------|----------|--------------------------------------|
| 007-77 | 200 µg | \$100.00 | <input type="button" value="ORDER"/> |

Product description

More information

[TLQP and other VGF-derived Peptides](#)

[Orphan Bioactive Peptides from VGF and PC2](#)

[Alzheimer biomarker](#)

| | |
|---------------------------|---|
| Catalog # | 007-77 |
| Standard Size | 200 µg |
| Sequence | Thr - Leu - Gln - Pro - Pro - Ala - Ser - Ser - Arg - Arg - Arg |
| Species | Rat |
| M.W | 1268.44 |
| Mass Spec Analysis | Exhibits correct M.W. |
| Solubility | Soluble in water. |
| Storage Condition | Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not refreeze any unused portions. |
| Appearance | White powder |
| Content | The contents of this vial have been accurately determined to contain 200µg of NET peptide. |
| HPLC Analysis | Purity ≥95% |