

TLQP-30 / VGF, Prepro (556-585) (Rat)

Catalog#	Standard size	Price	Order
007-79	100µg	\$250.00	<input type="button" value="ORDER"/>

Product description

More information

[TLQP and other VGF-derived Peptides](#)

[Orphan Bioactive Peptides from VGF and PC2](#)

[Alzheimer biomarker](#)

Catalog #	007-79
Standard Size	100µg
Sequence	Thr - Leu - Gln - Pro - Pro - Ala - Ser - Ser - Arg - Arg - Arg - His - Phe - His - His - Ala - Leu - Pro - Pro - Ala - Arg His - His - Pro - Asp - Leu - Glu - Ala Gln - Ala
Species	Rat
M.W	3745.58
Mass Spec Analysis	Exhibits correct M.W.
Storage Condition	Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not refreeze any unused portions.
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.

Thursday 20 December, 2012

11979368 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.