

TLQP-62 (Rat)

Catalog#	Standard size	Price	Order
007-82	100 µg	\$375.00	<div>ORDER</div>

- Product description
- More information

- TLQP and other VGF-derived Peptides
- Big NERP-2 & TLQP-62
- Orphan Bioactive Peptides from VGF and PC2
- Alzheimer biomarker

Catalog #	007-82
Standard Size	100 µg
Sequence	Thr - Leu - Gln - Pro - Pro - Ala - Ser - Ser - Arg - Arg - Arg - His - Phe - His - His - Ala - Leu - Pro - Pro - Ala - Arg - His - His - Pro - Asp - Leu - Glu - Ala - Gln - Ala - Arg - Arg - Ala - Gln - Glu - Glu - Ala - Asp - Ala - Glu - Glu - Arg - Arg - Leu - Gln - Glu - Gln - Glu - Glu - Leu - Glu - Asn - Tyr - Ile - Glu - His - Val - Leu - Leu - His - Arg - Pro
Species	Rat
M.W	7401.13
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.
Appearance	White Powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥95%