

## Prepro-Augurin (42-70) (human)

Catalog#	Standard size	Price	Order
012-19	100 µg	\$220.00	<input type="button" value="ORDER"/>

### Product description

### More information

[Augurin - Research Abstracts, Antibody Stainings, Kits Standard Curves, Sequences, etc.](#)

<b>Catalog #</b>	012-19
<b>Standard Size</b>	100 µg
<b>Sequence</b>	Glu - Ala - Pro - Val - Pro - Thr - Lys - Thr - Lys - Val - Ala - Val - Asp - Glu - Asn - Lys - Ala - Lys - Glu - Phe - Leu - Gly - Ser - Leu - Lys - Arg - Gln - Lys - Arg
<b>Species</b>	Human
<b>M.W</b>	3268.84
<b>Mass Spec Analysis</b>	Exhibits correct M.W.
<b>Solubility</b>	Soluble in Water
<b>Storage Condition</b>	Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.
<b>Appearance</b>	White Powder
<b>Content</b>	The contents of this vial have been accurately determined to contain 100ug of NET peptide
<b>HPLC Analysis</b>	> 95% Pure

Thursday 20 December, 2012

11979572 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.