

Prepro-Augurin (42-70) (mouse)

Catalog#	Standard size	Price	Order
012-20	100 µg	\$220.00	<input type="button" value="ORDER"/>

Product description

More information

[Augurin - Research Abstracts, Antibody Stainings, Kits Standard Curves, Sequences, etc.](#)

Catalog #	012-20
Standard Size	100 µg
Sequence	Glu - Gly - Pro - Val - Pro - Ser - Lys - Thr - Asn - Val - Ala - Val - Ala - Glu - Asn - Thr - Ala - Lys - Glu - Phe - Leu - Gly - Gly - Leu - Lys - Arg - Ala - Lys - Arg
Species	Mouse
M.W	3066.69
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water
Storage Condition	Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.
Appearance	White Powder
Content	The contents of this vial have been accurately determined to contain 100ug of NET peptide
HPLC Analysis	> 95% Pure

Thursday 20 December, 2012

11979574 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.