

Augurin (Human)

Catalog#	Standard size	Price	Order
012-25A	100 µg	\$380.00	<input type="button" value="ORDER"/>

Product description

More information

[Augurin - Research Abstracts, Antibody Stainings, Kits Standard Curves, Sequences, etc.](#)

[Isolation and Identification of Augurin Peptide from the Bovine Hypothalamus](#)

Catalog #	012-25A
Standard Size	100 µg
Sequence	pGlu - Leu - Trp - Asp - Arg - Thr - Arg - Pro - Glu - Val - Gln - Gln - Trp - Tyr - Gln - Gln - Phe - Leu - Tyr - Met - Gly - Phe - Asp - Glu - Ala - Lys - Phe - Glu - Asp - Asp - Ile - Thr - Tyr - Trp - Leu - Asn - Arg - Asp - Arg - Asn - Gly - His - Glu - Tyr - Tyr - Gly - Asp - Tyr - Tyr - Gln - Arg - His - Tyr - Asp - Glu - Asp - Ser - Ala - Ile - Gly - Pro - Arg - Ser - Pro - Tyr - Gly - Phe - Arg - His - Gly - Ala - Ser - Val - Asn - Tyr - Asp - Asp - Tyr
Species	Human
M.W	9690.43
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.
Appearance	White Powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥95%

Thursday 20 December, 2012

11979581 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.