

Augurin (Human) Short Form

Catalog#	Standard size	Price	Order
012-26A	100 µg	\$350.00	<input type="button" value="ORDER"/>

Product description

More information

[Augurin - Research Abstracts, Antibody Stainings, Kits Standard Curves, Sequences, etc.](#)

Catalog #	012-26A
Standard Size	100 µg
Sequence	pGlu - Leu - Trp - Asp - Arg - Thr - Arg - Pro - Glu - Val - Gln - Gln - Trp - Tyr - Gln - Gln - Phe - Leu - Tyr - Met - Gly - Phe - Asp - Glu - Ala - Lys - Phe - Glu - Asp - Asp - Ile - Thr - Tyr - Trp - Leu - Asn - Arg - Asp - Arg - Asn - Gly - His - Glu - Tyr - Tyr - Gly - Asp - Tyr - Tyr - Gln - Arg - His - Tyr - Asp - Glu - Asp - Ser - Ala - Ile - Gly - Pro - Arg
Species	Human
M.W	7843.64
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in buffer with pH 7.0-7.5.
Storage Condition	Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.
Appearance	White Powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥95%

Thursday 20 December, 2012

11979584 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.