

CRSP-1 (Porcine)

Catalog#	Standard size	Price	Order
015-62	100 µg	\$190.00	<input type="button" value="ORDER"/>

Product description

More information

References

Calcitonin Receptor-Stimulating Peptide -1, -2 & -3 (CRSP -1, -2 &-3)

CRSP-1

Catalog #	015-62
Standard Size	100 µg
Sequence	Ser - Cys - Asn - Thr - Ala - Thr - Cys - Met - Thr - His - Arg - Leu - Val - Gly - Leu - Leu - Ser - Arg - Ser - Gly - Ser - Met - Val - Arg - Ser - Asn - Leu - Leu - Pro - Thr - Lys - Met - Gly - Phe - Lys - Val - Phe - Gly - NH ₂
Species	Porcine
M.W	4099.35
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥96%