

Amylin (Rat, Mouse)

Catalog#	Standard size	Price	Order
017-11	500 µg	\$265.00	<input type="button" value="ORDER"/>

Product description

More information

References

Amylin

Normal Plasma Peptide/Protein Levels

Catalog #	017-11
Standard Size	500 µg
Sequence	Lys - Cys - Asn - Thr - Ala - Thr - Cys - Ala - Thr - Gln - Arg - Leu - Ala - Asn - Phe - Leu - Val - Arg - Ser - Ser - Asn - Asn - Leu - Gly - Pro - Val - Leu - Pro - Pro - Thr - Asn - Val - Gly - Ser - Asn - Thr - Tyr - NH ₂ Disulfide Bridge: Cys ₂ - Cys ₇
Species	Rat, Mouse
M.W	3920.45
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Insoluble in water. Rehydrate the entire vial with at least 100µl of 0.1M Acetic Acid or any other acidic buffer.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 500µg of NET peptide.
HPLC Analysis	Purity ≥95%
Additional articles	For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.