

Amyloid beta (11-42) (Human)

| Catalog# | Standard size | Price | Order |
|----------|---------------|----------|--------------------------------------|
| 018-22 | 100 µg | \$220.00 | <input type="button" value="ORDER"/> |

Product description

More information

Amyloid Beta Peptides

| | |
|--------------------|---|
| Catalog # | 018-22 |
| Standard Size | 100 µg |
| Sequence | Glu - Val - His - His - Gln - Lys - Leu - Val - Phe - Phe - Ala - Glu - Asp - Val - Gly - Ser - Asn - Lys - Gly - Ala - Ile - Ile - Gly - Leu - Met - Val - Gly - Gly - Val - Val - Ile - Ala |
| Species | Human |
| Mass Spec Analysis | Exhibits correct M.W. |
| Storage Condition | Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions. |
| Appearance | White powder |
| Content | The contents of this vial have been accurately determined to contain 100µg of NET peptide. |

Thursday 20 December, 2012

11979715 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.