

APL 1beta-28

Catalog#	Standard size	Price	Order
018-42	100 µg	\$200.00	<input type="button" value="ORDER"/>

Product description

More information

[Amyloid-beta Precursor-like Proteins \(APLPs\)](#)

[Peptide Fingerprint of Alzheimer s](#)

[Alzheimer's Disease fingerprint Peptide Library](#)

Catalog #	018-42
Standard Size	100 µg
Sequence	Asp - Glu - Leu - Ala - Pro - Ala - Gly - Thr - Gly - Val - Ser - Arg - Glu - Ala - Val - Ser - Gly - Leu - Leu - Ile - Met - Gly - Ala - Gly - Gly - Gly - Ser - Leu
M.W	2585.94
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100 µg of NET peptide.
HPLC Analysis	Purity ≥95%

Thursday 20 December, 2012

11979731 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.