

Quality Control Record

Catalog No.: 018-69

Lot No.: 427502

Product: APP Intracellular Domain (AICD) Short Form / β -Amyloid Precursor (727-770) (Human)

Sequence: pGlu-Tyr-Thr-Ser-Ile-His-His-Gly-Val-Val-Glu-Val-Asp-Ala-Ala-Val-Thr-Pro-Glu-Glu-Arg-His-Leu-Ser-Lys-Met-Gln-Gln-Asn-Gly-Tyr-Glu-Asn-Pro-Thr-Tyr-Lys-Phe-Phe-Glu-Gln-Met-Gln-Asn

Molecular Weight:
5165.67

Mass Spec Analysis:
Exhibits correct M.W.

HPLC Analysis:
Purity \geq 95%

Appearance: White powder

Peptide Content:
The contents of this vial have been accurately determined to contain 100 μ g of NET peptide.

Solubility and Reconstitution:
Soluble in water.

Storage Conditions:
Store peptide in DRY form at 0-5°C.

Notes: For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.

Quality Control Manager: _____ **Date:** _____
Galina Citowsky