

UN-3 (Human)

Catalog#	Standard size	Price	Order
019-43	100 µg	\$150.00	<input type="button" value="ORDER"/>

Product description

More information

References

[UN 1,2,3 / Thrombin Preproprotein and Comb 1](#)

Catalog #	019-43
Standard Size	100 µg
Sequence	Glu - Leu - Leu - Glu - Ser - Tyr - Ile - Asp - Gly - Arg - Pro - Thr - Ala - Thr - Ser - Glu - Tyr - Gln - Thr - Phe - Phe - Asn - Pro - Arg - NH ₂
Species	Human
M.W	2834.08
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in small amount of water.
Storage Condition	Store reconstituted polypeptide solution in 50 mM Tris-HCl, 10 mM reduced Glutathione, pH=8.0 in the elution buffer at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 10µg of NET peptide.
HPLC Analysis	Purity ≥95%
Additional articles	For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.

Thursday 20 December, 2012

11979790 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.