

## Comb-1 (Human)

Catalog#	Standard size	Price	Order
019-44	200 µg	\$80.00	<input type="button" value="ORDER"/>

### Product description

### More information

### References

[UN 1,2,3 / Thrombin Preproprotein and Comb 1](#)

<b>Catalog #</b>	019-44
<b>Standard Size</b>	200 µg
<b>Sequence</b>	Glu - Gly - Leu - Glu - Pro - Gly
<b>Species</b>	Human
<b>M.W</b>	600.63
<b>Mass Spec Analysis</b>	Exhibits correct M.W.
<b>Solubility</b>	Soluble in water.
<b>Storage Condition</b>	Store peptide in DRY form at 0-5°C.
<b>Appearance</b>	White powder
<b>Content</b>	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
<b>HPLC Analysis</b>	Purity ≥95%
<b>Additional articles</b>	For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.

Thursday 20 December, 2012

11979791 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.