

Galanin (1-16) (Rat, Porcine)

Catalog#	Standard size	Price	Order
026-08	200 µg	\$42.00	<input type="button" value="ORDER"/>

Product description

More information

References

Galanin and Its Antagonists

Neuron Function Modulation

Catalog #	026-08
Standard Size	200 µg
Sequence	Gly - Trp - Thr - Leu - Asn - Ser - Ala - Gly - Tyr - LeuLeu - Gly - Pro - His - Ala - Ile -
Species	Rat, Porcine
M.W	1668.87
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 200µg of NET peptide.
HPLC Analysis	Purity ;Ÿ97%