

GLP-2 (Human)

Catalog#	Standard size	Price	Order
028-15	100 µg	\$190.00	<input type="button" value="ORDER"/>

Product description

More information

References

GLP-1 & GLP-2

Normal Plasma Peptide/Protein Levels

Neuron Function Modulation

Catalog #	028-15
Standard Size	100 µg
Sequence	His - Ala - Asp - Gly - Ser - Phe - Ser - Asp - Glu - Met - Asn - Thr - Ile - Leu - Asp - Asn - Leu - Ala - Ala - Arg - Asp - Phe - Ile - Asn - Trp - Leu - Ile - Gln - Thr - Lys - Ile - Thr - Asp
Species	Human
M.W	3766.18
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Insoluble in water. Dissolve the peptide in 100 µl of 60% Acetonitrile in water with 0.1% TFA or DMSO, than dilute with any desired buffer.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥95%