

Glicentin-related Polypeptide (21-50) (Human)

Catalog#	Standard size	Price	Order
028-17	100 µg	\$190.00	<input type="button" value="ORDER"/>

Product description

References

Catalog #	028-17
Standard Size	100 µg
Sequence	Arg - Ser - Leu - Gln - Asp - Thr - Glu - Glu - Lys - Ser - Arg - Ser - Phe - Ser - Ala - Ser - Gln - Ala - Asp - Pro - Leu - Ser - Asp - Pro - Asp - Gln - Met - Asn - Glu - A
Species	Human
M.W	3384.53
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥95%

Thursday 20 December, 2012

11980348 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.