

JI-38

//-->

Product description

More information

JI-38 - Potent GHRH agonist

Catalog #	031-03
Standard Size	100 µg
Species	Human
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥95%

Tuesday 29 March, 2011

6936853 requests since Wednesday 05 April, 2006

Copyright © 2011 PHOENIX PHARMACEUTICALS, INC.