

## Dap3 (Octanoyl)-Ghrelin (Human)

Catalog#	Standard size	Price	Order
031-58	100 µg	\$220.00	<input type="button" value="ORDER"/>

### Product description

### More information

### References

[All About Ghrelin](#)

[Ghrelin Kits - more information](#)

[Normal Plasma Peptide/Protein Levels](#)

<b>Catalog #</b>	031-58
<b>Standard Size</b>	100 µg
<b>Sequence</b>	Gly - Ser - Dap3(Octanoyl) - Phe - Leu - Ser - Pro - Glu - His - Gln - Arg - Val - Gln - Gln - Arg - Lys - Glu - Ser - Lys - Lys - Pro - Pro - Ala - Lys - Leu - Gln - Pro - Arg
<b>Species</b>	Human
<b>M.W</b>	3369.92
<b>Mass Spec Analysis</b>	Exhibits correct M.W.
<b>Solubility</b>	Soluble in water.
<b>Storage Condition</b>	Store peptide in DRY form at 0-5°C.
<b>Appearance</b>	White powder
<b>Content</b>	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
<b>HPLC Analysis</b>	Purity ≥95%
<b>Additional articles</b>	more stable than active Ghrelin in solution - same biological finction - web shows grapfic

Thursday 20 December, 2012

11980522 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.