

Ghrelin (1-5) Amide, [Dap3]-Octanoyl (Human, Rat, Mouse)

| Catalog# | Standard size | Price | Order |
|----------|---------------|----------|--------------------------------------|
| 032-14 | 200 µg | \$130.00 | <input type="button" value="ORDER"/> |

Product description

More information

References

[All About Ghrelin](#)

[Ghrelin Kits - more information](#)

[Normal Plasma Peptide/Protein Levels](#)

| | |
|----------------------------|--|
| Catalog # | 032-14 |
| Standard Size | 200 µg |
| Sequence | Gly - Ser - Dap3(Octanoyl) - Phe - Leu - NH ₂ |
| Species | Human, Rat, Mouse |
| M.W | 631.72 |
| Mass Spec Analysis | Exhibits correct M.W. |
| Solubility | Soluble in water. |
| Storage Condition | Store peptide in DRY form at 0-5°C. |
| Appearance | White powder |
| Content | The contents of this vial have been accurately determined to contain 200 µg of NET peptide. Do not attempt to weigh out smaller portions. |
| HPLC Analysis | Purity ≥96% |
| Additional articles | For best results and reproducibility, rehydrate peptide just before use. Do not refreeze any unused portions. Do not attempt to weigh out smaller portions |

Thursday 20 December, 2012

11980576 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.