

GMA(Dap3-Oct)-Ghrelin-NH2 (Rat, Mouse)

Catalog#	Standard size	Price	Order
032-34	200 µg	\$220.00	<input type="button" value="ORDER"/>

Product description

More information

[All About Ghrelin](#)

[Ghrelin Kits - more information](#)

Normal Plasma Peptide/Protein Levels

Catalog #	032-34
Standard Size	200 µg
Sequence	Gly - Met - Ala - Dap3(Octanyl) - Phe - Leu - Ser - Pro - Glu - His - Gln - Lys - Ala - Gln - Gln - Arg - Lys - Glu - Ser - Lys - Lys - Pro - Pro - Ala - Lys - Leu - Gln - Pro - Arg - NH2
Species	Rat, Mouse
M.W	3573.19
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 200µg of NET peptide.
HPLC Analysis	Purity ≥95%

Thursday 20 December, 2012

11980578 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.