

INSL 3 C-Peptide (Human)

Catalog#	Standard size	Price	Order
035-48	100 µg		<input type="button" value="ORDER"/>

Product description

More information

References

[Relaxin & Relaxin-Related Peptides](#)

[C-Peptide from Proinsulin](#)

[Insulin-like 3 Peptide \(INSL 3\)](#)

[C-Peptide from Prepro-Insulin-like Peptide](#)

[Relaxin Evolution](#)

Catalog #	035-48
Standard Size	100 µg
Sequence	Pro - Ala - Ala - Gly - Gly - Asp - Arg - Glu - Leu - Leu - Gln - Trp - Leu - Glu - Arg - Arg - His - Leu - Leu - His - Glu - Leu - Val - Ala - Asp - Ser - Asn - Leu - Thr - Leu - Gly - Pro - Gly - Leu - Gln - Pro - Leu - Pro - Gln - Thr - Ser - His - His - His - Arg - His - His - Arg
Species	Human
M.W	5426.12
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥ 95%

Thursday 20 December, 2012

11980677 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.