

P-94-Gly Amide (Human)

Catalog#	Standard size	Price	Order
036-35	100 µg	\$350.00	<input type="button" value="ORDER"/>

- Product description
- More information
- References

New Peptide Agonists for GPCRs

Catalog #	036-35
Standard Size	100 µg
Sequence	Ala - Ala - Gln - Ala - Thr - Gly - Pro - Leu - Gln - Asp - Asn - Glu - Leu - Pro - Gly - Leu - Asp - Glu - Arg - Pro - Pro - Arg - Ala - His - Ala - Gln - His - Phe - His - Lys - His - Gln - Leu - Trp - Pro - Ser - Pro - Phe - Arg - Ala - Leu - Lys - Pro - Arg - Gly - NH2
Species	Human
M.W	5170.84
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100 µg of NET peptide
HPLC Analysis	Purity ≥95%