

LH-RH / Gn-RH (Human, Rat, Mouse, Porcine)

Catalog#	Standard size	Price	Order
040-02	500 µg	\$27.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Prepro-Luteinizing Hormone-Releasing Hormone \(Human\) Prohormone Schematic](#)

Catalog #	040-02
Standard Size	500 µg
Sequence	pGlu - His - Trp - Ser - Tyr - Gly - Leu - Arg - Pro - Gly - NH ₂
Species	Human, Rat, Porcine, Mouse
M.W	1182.31
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 500µg of NET peptide.
HPLC Analysis	Purity ≥96%
Additional articles	For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.

Thursday 20 December, 2012

11980767 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.