

Neuromedin U-9 (Mouse, Guinea Pig)

Catalog#	Standard size	Price	Order
046-37	200 µg	\$130.00	<input type="button" value="ORDER"/>

Product description

More information

References

Neuromedin U (NmU)

Catalog #	046-37
Standard Size	200 µg
Sequence	Gly - Tyr - Phe - Leu - Phe - Arg - Pro - Arg - Asn - NH ₂
Species	Mouse
Specific Specie	Guinea Pig
M.W	1168.37
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 200µg of NET peptide.
HPLC Analysis	Purity ≥96%
Additional articles	For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.