

Neuromedin U-8 (Porcine)

| Catalog# | Standard size | Price | Order |
|----------|---------------|---------|--------------------------------------|
| 046-39 | 500 µg | \$32.00 | <input type="button" value="ORDER"/> |

Product description

More information

References

Neuromedin U (NmU)

| | |
|----------------------------|---|
| Catalog # | 046-39 |
| Standard Size | 500 µg |
| Sequence | Tyr - Phe - Leu - Phe - Arg - Pro - Arg - Asn - NH ₂ |
| Species | Porcine |
| M.W | 1110.59 |
| Mass Spec Analysis | Exhibits correct M.W. |
| Solubility | Soluble in water. |
| Storage Condition | Store peptide in DRY form at 0-5°C. |
| Appearance | White powder |
| Content | The contents of this vial have been accurately determined to contain 500µg of NET peptide. |
| HPLC Analysis | Purity ≥95% |
| Additional articles | For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions. |