

## Neuromedin U-25 (Porcine)

Catalog#	Standard size	Price	Order
046-40	200 µg	\$42.00	<input type="button" value="ORDER"/>

### Product description

### More information

### References

#### Neuromedin U (NmU)

<b>Catalog #</b>	046-40
<b>Standard Size</b>	200 µg
<b>Sequence</b>	Phe - Lys - Val - Asp - Glu - Glu - Phe - Gln - Gly - Pro - Ile - Val - Ser - Gln - Asn - Arg - Arg - Tyr - Phe - Leu - Phe - Arg - Pro - Arg - Asn - NH <sub>2</sub>
<b>Species</b>	Porcine
<b>M.W</b>	3140.63
<b>Mass Spec Analysis</b>	Exhibits correct M.W.
<b>Solubility</b>	Soluble in water.
<b>Storage Condition</b>	Store peptide in DRY form at 0-5°C.
<b>Appearance</b>	White powder
<b>Content</b>	The contents of this vial have been accurately determined to contain 200µg of NET peptide.
<b>HPLC Analysis</b>	Purity ≥96%

Thursday 20 December, 2012

11980962 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.