

Neuromedin U-25 (Human)

Catalog#	Standard size	Price	Order
046-42	200 µg	\$160.00	<input type="button" value="ORDER"/>

Product description

More information

References

Neuromedin U (NmU)

Catalog #	046-42
Standard Size	200 µg
Sequence	Phe - Arg - Val - Asp - Glu - Glu - Phe - Gln - Ser - Pro - Phe - Ala - Ser - Gln - Ser - Arg - Gly - Tyr - Phe - Leu - Phe - Arg - Pro - Arg - Asn - NH ₂
Species	Human
M.W	3080.42
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 200µg of NET peptide.
HPLC Analysis	Purity ≥96%
Additional articles	For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.