

Neuromedin U, Prepro (106-141) (Rat, Mouse)

| Catalog# | Standard size | Price | Order |
|----------|---------------|----------|--------------------------------------|
| 046-54 | 100 µg | \$190.00 | <input type="button" value="ORDER"/> |

Product description

More information

References

Neuromedin U (NmU)

| | |
|---------------------------|---|
| Catalog # | 046-54 |
| Standard Size | 100 µg |
| Sequence | Phe - Leu - Phe - His - Tyr - Ser - Lys - Thr - Gln - Lys - Leu - Gly - Asn - Ser - Asn - Val - Val - Ser - Ser - Val - Val - His - Pro - Leu - Leu - Gln - Leu - Val - Pro - Gln - Leu - His - Glu - Arg - Arg - Met |
| Species | Rat, Mouse |
| M.W | 4203.93 |
| Mass Spec Analysis | Exhibits correct M.W. |
| Solubility | Soluble in water. |
| Storage Condition | Store peptide in DRY form at 0-5°C. |
| Appearance | White powder |
| Content | The contents of this vial have been accurately determined to contain 100µg of NET peptide. |
| HPLC Analysis | Purity ≥95% |

Thursday 20 December, 2012

11980975 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.