

## Neuromedin U-23 (Mouse)

Catalog#	Standard size	Price	Order
046-65	100 µg	\$190.00	<input type="button" value="ORDER"/>

### Product description

### More information

### References

#### Neuromedin U (NmU)

<b>Catalog #</b>	046-65
<b>Standard Size</b>	100 µg
<b>Sequence</b>	Phe - Lys - Ala - Glu - Tyr - Gln - Ser - Pro - Ser - Val - Gly - Gln - Ser - Lys - Gly - Tyr - Phe - Leu - Phe - Arg - Pro - Arg - Asn - NH <sub>2</sub>
<b>Species</b>	Mouse
<b>M.W</b>	2706.06
<b>Mass Spec Analysis</b>	Exhibits correct M.W.
<b>Solubility</b>	Soluble in water.
<b>Storage Condition</b>	Store peptide in DRY form at 0-5°C.
<b>Appearance</b>	White powder
<b>Content</b>	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
<b>HPLC Analysis</b>	Purity ≥95%
<b>Additional articles</b>	Howard et al. Identification of receptors for neuromedin U and its role in feeding. Nature. 2000 Jul 6;406(6791):70-4. For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.