

Neuromedin U, Prepro (104-139) (Human)

Catalog#	Standard size	Price	Order
046-76	100 µg	\$160.00	<input type="button" value="ORDER"/>

Product description

More information

References

Neuromedin U (NmU)

Catalog #	046-76
Standard Size	100 µg
Sequence	Phe - Leu - Phe - His - Tyr - Ser - Lys - Thr - Gln - Lys - Leu - Gly - Lys - Ser - Asn - Val - Val - Ser - Ser - Val - Val - His - Pro - Leu - Leu - Gln - Leu - Val - Pro - His - Leu - His - Glu - Arg - Arg - Met
Species	Human
M.W	4227.01
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity =95%

Thursday 20 December, 2012

11981018 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.