

QRFP (28n) (Human)

Catalog#	Standard size	Price	Order
049-59	100 µg	\$160.00	<input type="button" value="ORDER"/>

Product description

More information

References

QRFP-43

Catalog #	049-59
Standard Size	100 µg
Sequence	pGlu - Asp - Glu - Gly - Ser - Glu - Ala - Thr - Gly - Phe - Leu - Pro - Ala - Ala - Gly - Glu - Lys - Thr - Ser - Gly - Pro - Leu - Gly - Asn - Leu - Ala - Glu - Glu
Species	Human
M.W	2759.91
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of peptide.
HPLC Analysis	Purity ≥95%

Thursday 20 December, 2012

11981360 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.