

QRFP 43 (Mouse)

Catalog#	Standard size	Price	Order
049-61	100 µg	\$235.00	<input type="button" value="ORDER"/>

Product description

More information

References

QRFP-43

Catalog #	049-61
Standard Size	100 µg
Sequence	pGlu - Asp - Gly - Ser - Ser - Glu - Ala - Ala - Gly - Phe - Leu - Pro - Ala - Asp - Ser - Glu - Lys - Ala - Ser - Gly - Pro - Leu - Gly - Thr - Leu - Ala - Glu - Glu - Leu - Ser - Ser - Tyr - Ser - Arg - Arg - Lys - Gly - Gly - Phe - Ser - Phe - Arg - Phe - NH ₂
Species	Mouse
M.W	4493.90
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Insoluble in water. Dissolve peptide in 50 µl of 60% Acetonitrile in water with 0.1% TFA or DMSO, then dilute with any desired buffer.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of peptide.
HPLC Analysis	Purity ≥95%

Thursday 20 December, 2012

11981362 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.