

BRAWNIN Delta-25 / C12orf73 (26-71) (Human)

Quantity:

Order

Catalog #	Size	Price
051-25	100 µg	\$250

Technical

Summary

Sequence: (Ala⁵A)

Glu-Val-Val-His-Arg-Tyr-Tyr-Arg-Pro-Asp-Leu-Thr-
 Ile-Pro-Glu-Ile-Pro-Pro-Lys-Arg-Gly-Glu-Leu-Lys-
 Thr-Glu-Leu-Leu-Gly-Leu-Lys-Glu-Arg-Lys-His-Lys-
 Pro-Gln-Val-Ser-Gln-Gln-Glu-Glu-Leu-Lys

Molecular Weight:	5507.31
Purity:	≥ 95%
Validation:	Exhibits correct molecular weight
Solubility:	Soluble in water
Storage:	Up to 6 months in lyophilized form at 0-5°C. For best results, rehydrate just before use. Aliquot before freezing to avoid repeated freeze-thaw cycles.
Appearance:	White powder
Contents:	Each vial contains 100 µg of NET peptide.



(<https://phoenixpeptide.com/cemark>)



(<https://phoenixpeptide.com/iso>)

About Us (<https://phoenixpeptide.com/about>) Privacy Policy (<https://phoenixpeptide.com/privacy>)

Terms of Use (<https://phoenixpeptide.com/terms>) (https://phoenixpeptide.com/order_information)

Copyright © 2016 Phoenix Pharmaceuticals.

ABOUT SSL CERTIFICATES

(<https://www.websecurity.symantec.com/ssl-certificate>)