

Exendin-3 (9-39)-Amide (Heloderma Horridum)

Catalog#	Standard size	Price	Order
070-93	100 µg	\$105.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Exendin-4, Research Abstracts, Kits Standard Curves, Sequences, etc.](#)

[Neuron Function Modulation](#)

Catalog #	070-93
Standard Size	100 µg
Sequence	Asp - Leu - Ser - Lys - Gln - Met - Glu - Glu - Glu - Ala - Val - Arg - Leu - Phe - Ile - Glu - Trp - Leu - Lys - Asn - Gly - Gly - Pro - Ser - Ser - Gly - Ala - Pro - Pro - Pro - Ser - NH ₂
Specific Specie	Heloderma horridum
M.W	3367.67
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥95%
Additional articles	For best results and reproducibility, rehydrate peptide just before use. Do not refreeze any unused portions. Do not attempt to weigh out smaller portions.

Thursday 20 December, 2012

11982146 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.