

Exendin-4 (Heloderma suspectum)

Catalog#	Standard size	Price	Order
070-94	100 µg	\$140.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Exendin-4](#), [Research Abstracts](#), [Kits](#) [Standard Curves](#), [Sequences](#), etc.

Neuron Function Modulation

Catalog #	070-94
Standard Size	100 µg
Sequence	His - Gly - Glu - Gly - Thr - Phe - Thr - Ser - Asp - Leu - Ser - Lys - Gln - Met - Glu - Glu - Glu - Ala - Val - Arg - Leu - Phe - Ile - Glu - Trp - Leu - Lys - Asn - Gly - Gly - Pro - Ser - Ser - Gly - Ala - Pro - Pro - Pro - Ser - NH ₂
Specific Specie	Heloderma suspectum
M.W	4184.01
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥95%
Additional articles	For best results and reproducibility, rehydrate peptide just before use. Do not refreeze any unused portions. Do not attempt to weigh out smaller portions.