

Exendin-4 (3-39) (Heloderma suspectum) / P1, Inhibitor of Platelet Calpain

Catalog#	Standard size	Price	Order
070-99	100 µg	\$130.00	<input type="button" value="ORDER"/>

Product description

More information

References

[Exendin-4, Research Abstracts, Kits Standard Curves, Sequences, etc.](#)

Neuron Function Modulation

Catalog #	070-99
Standard Size	100 µg
Sequence	Glu - Gly - Thr - Phe - Thr - Ser - Asp - Leu - Ser - Lys - Gln - Met - Glu - Glu - Glu - Ala Val - Arg - Leu - Phe - Ile - Glu - Trp - Leu - Lys - Asn - Gly - Gly - Pro - Ser - Ser - Gly - Ala - Pro - Pro - Pro - Ser - NH2
Specific Specie	Heloderma suspectum
M.W	3992.44
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥95%

Thursday 20 December, 2012

11982152 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.