

LEAP-2 (38-77) (Human)

Catalog#	Standard size	Price	Order
075-40	100 µg	\$235.00	<input type="button" value="ORDER"/>

Product description

More information

References

Anti-Microbial Peptides

LEAP-2

Catalog #	075-40
Standard Size	100 µg
Sequence	Met - Thr - Pro - Phe - Trp - Arg - Gly - Val - Ser - Leu - Arg - Pro - Ile - Gly - Ala - Ser - Cys - Arg - Asp - Asp - Ser - Glu - Cys - Ile - Thr - Arg - Leu - Cys - Arg - Lys - Arg - Arg - Cys - Ser - Leu - Ser - Val - Ala - Gln - Glu ~ Disulfide bridges between Cys54 & Cys 65 and Cys60 & Cys70
Species	Human
M.W	4583.35
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity ≥95%

Thursday 20 December, 2012

11982300 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.