

Phoenixin-20 amide (Human, Rat, Mouse, Bovine)

Catalog#	Standard size	Price	Order
079-03	200ug	\$150.00	<input type="button" value="ORDER"/>

Product description

More information

Phoenixin - Exclusively from Phoenix Pharmaceuticals

Catalog #	079-03
Standard Size	200ug
Sequence	Ala - Gly - Ile - Val - Gln - Glu - Asp - Val - Gln - Pro - Pro - Gly - Leu - Lys - Val - Trp - Ser - Asp - Pro - Phe - NH ₂
Species	Human, Rat, Mouse, Bovine
M.W	2181.5
Mass Spec Analysis	Exhibit correct M.W.
Solubility	Soluble in water
Storage Condition	Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions
Appearance	White Powder
Content	200ug NET peptide
HPLC Analysis	> 95% purity

Thursday 20 December, 2012

11982528 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.