

OBAP (Human) - Biotin Labeled

Catalog#	Standard size	Price	Order
B-055-52	100 ug	\$350.00	<input type="button" value="ORDER"/>

Product description

More information

Osteoblast Activating Peptide (OBAP)

Catalog #	B-055-52
Standard Size	100 ug
Sequence	Biotinyl - Leu - Asp - Leu - Asn - Leu - Glu - Leu - Ser - Lys - Phe - Arg - Met - Pro - Gln - Pro - Ser - Ser - Gly - Arg - Glu - Ser - Pro - Arg - His
Species	Human
M.W	3021.68
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Before using, store the peptide in the DRY form at 0-5°C. For best and repeatable result, rehydrate peptide right before using. Do not re-freeze any unused portion.
Content	This vial contains 100µg of NET peptide.
Counter_Ion	TFA salt
HPLC Analysis	Purity ≥ 95%

Thursday 20 December, 2012

11982920 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.